

## Electrical Delivery Bed



- Fully electrical delivery bed.
- With easy detachable PP head and foot boards.
- PP side rails on head side of bed, angle indicator on each side
- The foot section is easily removed and installed,

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Fully electrical delivery bed.

- With easy detachable PP head and foot boards.
  - PP side rails on head side of bed, angle indicator on each side
  - The foot section is easily removed and installed, serum hanger housing at - Four Corners and 12 cm diameter impact bumper wheels. Foot supports on
  - both sides steel tube chrome plated and straining handle.
  - High density foam mattress, liquid proof coated with vinlex, easy to - disinfectant clean.
  - Delivery bed backrest, height and Trendelenburg movements can be - adjusted easily and quietly with hand control.
  - Ability to move without difficulty and quietly by motor
  - Bed is mobile and can be moved easily through 4 rotatable wheels, two of - the wheels are braked and the brake is activated when necessary.
  - There is a stainless waste bowl in the delivery bed.
- Included with the bed is a height adjustable IV pole.

## BEDS & STRETCHERS: Electrical Delivery Bed

---

- The knee supports are with height adjustment chrome plated steel and the - footrest with upholstery top waterproof vinlex coating
- Bed safe working load is 200 kg
- Specifications - Back rest and Trendelenburg, anti Trendelenburg movements make by electrical motor
- Detachable PP mono block moulded head & foot boards
- Clean and smooth appearance
- Wall protective soft four corner bumpers Four Brackets on corner
- Strong and durable structure made of steel profile
- 360° rotatable 125mm diameter wheels two of them with brakes
- I.V pole with height adjustment