

## Bariatric Electric 5 Function Bed



- Headrest and height adjustment with electric motors.
- Fowler and vascular position.
- Hand control unit.
- Knee adjustment with gradual mechanism.
- Head and foot boards made from wood.
- 2 Pcs metal side rails. Metal bed surface.
- 360° movable castors with a linear braking system.
- Battery Backup.
- Overbed lifting pole.
- Low height. 290mm – 665mm Weight capacity 350 kg
- Mattress width 1090 mm
- Bed total width 1200 mm
- Bed length 2130 mm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)