

## Hydraulic Delivery Obstetrics Delivery Bed



Size: 2035/1510mm(L)\*940/664mm(W)\*665/950mm(H)

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description Main feature:

- Epoxy coated steel bed frame, with ABS engineering plastic side and head board;
- Telescopic working table, to be controlled by linear bearing, telescopic length : 550mm.
- Three functions as follow:
  - 1) Back rest:0-70±5° operate by gas spring.
  - 2) Height adjustable from 665-950mm, operate by hydraulic pedal.
  - 3) Trendelenburg & reverse Trendelenburg 14±2°, operate by gas spring.
- Equipped with high density sponge mattress covered with artificial leather, and the body part - thickness is 8cm and the telescopic part thickness is 14cm.
- Adjustable leg support and sliding grip holder.
- Four dia-125mm castors