

## Dual CPR Training Manikin



- Uniquely designed dual-function Adult and Child CPR Training Torso

Rating: Not Rated Yet  
**Price**

[Ask a question about this product](#)

Description - Switch between modes with manual dial

- Dial adjusts both the depth of compressions and volume of ventilations to that of either an Adult or a Child
- Saves you the need of having to buy separate manikins to teach both Adult and Child CPR
- Can be used with an AED Trainer
- One-way valve system on airway promotes excellent hygiene
- Replaceable face and chest skin
- Extra skins can be purchased as an addition to your Practiman and can be interchanged while cleaning
- Excellent quality
- Supplied with spare lungs and airway valves, as well as CPR Face Shields and Bag