

## Hydraulic Hoist



- Base on 4 wheels allows the hoist to be moved and used in different rooms
- Brakes on 2 wheels easy to apply and release
- Removable mast and boom for disassembly and transport
- Base spreading adjustment handle spreads the base legs for maximum stability- 5 settings

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description Supplied with a standard sling and spreader bar.

The mechanism is operated with a hydraulic system designed to assist with heavy lifting.

Safely lifts to a maximum user weight of 150kg.

Wide variety of slings are available to allow maximum comfort for the user.