



## Suitcase Ramp



- One piece aluminium
- Carrying capacity of 270kg
- Weight – 13.5kg
- Centre fold design
- Max Length – 1520mm
- Max step height – 15.2cm
- Sturdy handle
- 2 sets of fixation bolts & chains

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description This Suitcase Ramp is a one-piece aluminium ramp which can be used for a range of assistive devices. Its 1520mm length will allow you to access buildings with steps, kerbs and other obstacles.

The non-slip surface on the underside of the ramp prevents it from moving off the step or kerb.

The ramp folds in half and can be bolted shut to prevent it from opening.

A handle allows easy and convenient transport of the ramp.

These bolts are fixed with a small chain to ensure that they are not lost while the hoist is in use.

The ramp is also supplied with 2 sets of fixing bolts- Allowing you to secure the ramp permanently.