

## Telescopic Ramp



- Two piece aluminium
- 1.07m or 1.83m in length
- Each ramp is 200mm wide
- Product Weight – 10kg
- Carrying Capacity – 200kgs

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description These aluminium ramps have a weight capacity of 200kg and yet only weigh 10kg. They can be compressed and packed away in a small, convenient bag. The ramps can be extended from 1070mm to 1830mm, allowing you to get wheelchairs up and over stairs, curbs and other obstacles. There is an anti-slip coating on the top lip surface of the ramp, preventing it from sliding while being used