

Vacuum Mattress & Accessories



The Vacuum Mattress is the golden standard for full body immobilisation during transportation

- Manufactured from high quality TPU material
- 8 built in handles to ensure safe and stable handling of patient
- 4 built in patient restraining straps
- Smooth surface promotes hygiene and infection control during cleaning
- Vacuum Mattress is X-ray and CT scan friendly
- Double-chamber pump and bag included
- Particularly valuable in situations where:
 - There is a high degree of suspicion for spinal/vertebral injury
 - Patient has pelvic or extremity fractures
 - Immobilisation is required to minimise movement in the unstable intubated patient during transportation
- Load-bearing: 159kg
- Dimensions: • 92 × 53 × 140cm • 6kg
- Standards Compliance: • CE, FDA & ISO13485 approved

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

