


## Mattress Modular Cut



Modular LASER cut foam mattress with waterproof cover  
Size 1980 x 880 x 160 mm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description Modular LASER cut foam mattress with waterproof cover  
Size 1980 x 880 x 160 mm

Mattresses included here has a cross cut upper surface.

This results in the mattress somewhat conforming to the shape and movement of the body. This can reduce friction and shear by allowing the surface to move with the user.

A cross cut surface can also provide better ventilation for the skin as flat surfaced foam can act as an insulator and increase