

## Advanced Tube Mattress



### Tube Mattress with pump

- 22 Separate Tubes
- Navy Blue
- 200x90x10mm
- Cycle time: 12 minutes
- Supplied with Pump
- End flap included
- Waterproof cover included
- Supplied with a low pressure, high flow pump
- Tubes can be replaced if damaged (no need to throw the entire mattress away)
- Pressure-selector dial to mattress pressure
- End flaps to secure to bed mattress
- Optional-use waterproof cover
- Max user weight of 135kg

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)