

## Cruiser Wheelchair



• Size : 14" to 20" • Frame : Steel / Folding • Upholstery : Vinyl • Armrest : Removable / Full • Footrest : Swing-Away • Footplate : ABS • Castor : 8" x 2" Solid PU or 8" x 1" • Rear Wheel : 24" x 1.3/8" Mag Solid PU • Max. Carry Capacity : 110kgs • Weight of Chair : 18.8kgs • Back Height : 42cm • Seat Depth : 40cm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)