

## 5 Function Electrical Bed



### 5 Function Electrical Homecare Bed

- › Electrically operated backrest, height and knee break angle - Comfort sitting position for sitting upright
- › Four-part reclining surface with back rest adjustment
- › Infinitely variable, electrically operated head and foot tilt functionality (Trendelenburg & reverse Trendelenburg)
- › Auto contour-simultaneous adjustment of the backrest and knee break
  
- › Patient Lifting Pole optional
- › Integral full length steel or wooden side rails
- › Double Head and Foot Boards
- › Dewert Actuators
- › Silent Double Castors with brake
- › Integrated bed extension available
- › Extensive customisation options available
- › Easy to assemble, disassemble and transport

Rating: Not Rated Yet

**Price**

## BEDS: 5 Function Electrical Bed

---

[Ask a question about this product](#)

Description	MSCB1001
Overall length	216 cm
Overall width	100 cm
Height adjustment	32 to 71 cm
Mattress Platform length	202 cm
Mattress Platform width	87 cm
Backrest Tilt	70 degrees
Footrest Tilt	26 degrees
Safe working load	213 kg
Maximum user weight	180 kg