

Infant Cot Bed Bassinet



- Luxurious infant Bed (with Trendelenburg and reverse Trendelenburg 12°)
- Size: 870×530×920~1070mm

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description Main feature:

- Clear plastic bassinet to secure clear visibility for observer.
- Inclination angle adjustable by gas spring, 0-12°. Height adjustable by gas spring, 920-1070mm.
- 100mm dia. Polyurethane single-wheel casters with brakes on opposite corner, The smooth rotating caster will provide silent transportation. - blue colour available