

## Invacare Lite Wheelchair



- Size : 16" & 18"
- Frame : Blue - Lightweight / Folding

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description • Upholstery : Nylon

- Armrest : Removable / Full
- Footrest : Swing-Away
- Footplate : ABS
- Castor : 8" x 2" Solid PU
- Rear Wheel : 12.5"
- Max. Carry Capacity : 100kgs
- Weight of Chair : 12kgs
- Back Height : 45cm, Folding backrest
- Seat Depth : 43cm