

## Lightweight Breezy Wheelchair



- Size : 16?, 18? & 20?
- Frame : Aluminium /Folding

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description • Upholstery : Vinyl

- Armrest : Flip Back Desk/Height Adj. PU
- Footrest : Removable
- Rear Wheel : 24? x 1.3/8?
- Quick Release Back Wheel
- Max. Carry Capacity : 125Kgs
- Weight of Chair : 13.5kgs
- 5 Rear wheel axle positions for adjustments
- Adjustable Seat height and depth