

## Lightweight Half Back Wheelchair



• Size : 16", 18" & 20" • Frame : Aluminium Black/Folding

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description • Upholstery : Nylon • Armrest : Flip Back Desk/Height Adj. PU • Footrest : Removable/Aluminium Hanger • Footplate : ABS • Castor : 7" x 2" Solid PU • Rear Wheel : 24" x 1.3/8" • Quick Release Back Wheel • Half Back • Max. Carry Capacity : 100Kgs • Weight of Chair : 16kgs