

## Lightweight Wheelchair



- Size : 14" - 20"
- Frame : Lightweight / Folding

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description • Upholstery : Vinyl

- Armrest : Removable / Full
- Footrest : Swing-Away
- Footplate : ABS
- Castor : 8" x 2" Solid PU
  - Rear Wheel : Quick Release 24" x 1 3/8 Spoke Solid PU
- Max. Carry Capacity : 95kgs
- Weight of Chair : 14.2kgs
- Back Height : 44cm
- Seat Depth : 42cm