

## Manual Steel Wheelchair



• Size : 14" , 18" & 20" • Frame : Steel / Folding • Upholstery : Nylon • Armrest : Removable / Desk • Footrest : Removable / Steel hanger • Footplate : ABS • Castor : 8" x 2" Solid PU • Rear wheel : 24" x 1.3/8" Spoke Solid PU • Max. Carry Capacity : 80kgs – 90kgs • Weight of Chair : 21kgs

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)