

## Infant Cot Bed Bassinet



- Luxurious infant Bed (with Trendelenburg and reverse Trendelenburg 12°)
- Size: 870x530x920-1070mm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description Main feature:

- Clear plastic bassinet to secure clear visibility for observer.
- Inclination angle adjustable by gas spring, 0-12°. Height adjustable by gas spring, 920-1070mm.
- 100mm dia. Polyurethane single-wheel casters with brakes on opposite corner, The smooth rotating caster will provide silent transportaion. - blue colour available