

## Tuffee Wheelchair



- Size : 14" to 20"
- Frame : Steel / Folding
  - Upholstery : Reinforced Vinyl
- Armrest : Removable / Full
- Footrest : Swing-Away
- Footplate : ABS
- Castor : 8" x 2" Solid PU
- Rear Wheel : 24" x 1.75 Mag Solid PU
- Max. Carry Capacity : 110kgs
- Weight of Chair : 19.3kgs
- Back Height : 42cm
- Seat Depth : 41cm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)