

## Rainbow Wheelchair



- Size : 18"
- Frame : Blue / Steel / Folding
  - Upholstery : Vinyl
- Armrest : Removable / Full
- Footrest : Swing-Away
- Footplate : ABS
- Castor : 8" x 1" Solid PU
  - Rear Wheel : Fixed 24" x 1 3/8 Mag Solid PU
- Max. Carry Capacity : 100kgs
- Weight of Chair : 18.8kgs
- Back Height : 42cm
- Seat Depth : 40cm

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)