

Trauma Stretcher



Rating: Not Rated Yet

Price

[Ask a question about this product](#)

- Description - Back rest: $0-85^{\circ}\pm 5^{\circ}$, controlled by gas spring
- Lowest height: $600\text{mm}\pm 10\text{mm}$
 - Highest height: $920\text{mm}\pm 10\text{mm}$
 - Trendelenburg tilting back & forth is controlled by the foot pedal on the sides of the stretcher.
 - Safety working load: 2500N
 - Using environment : temperature: $5^{\circ}\text{--}40^{\circ}$
 - The main frame of stretcher is made of galvanized steel tube EPOXY coated electrostatic spray.
 - The top board BACKREST is made of X-ray material.