

## Trauma Stretcher

PRODUCT  
Image  
Unavailable



Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Back rest:  $0-85^{\circ}\pm 5^{\circ}$ , controlled by gas spring.

- Knee break:  $0-30^{\circ}\pm 5^{\circ}$ , controlled by manual crank

- Lowest height:  $600\text{mm}\pm 10\text{mm}$

- Highest height:  $920\text{ mm}\pm 10\text{mm}$

- Tilt to back & forth which is controlled by the foot pedal on the two sides of the stretcher.

- Safety working loading: 2500N

- User environment : temperature:  $5^{\circ}$ -  $40^{\circ}$

- The main frame of stretcher is made of galvanized steel tube coated with electrostatic powder spray.

- The top board is made of X-ray translucent material.