

Trauma Stretcher



Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description - Back rest: $0-85^{\circ} \pm 5^{\circ}$, controlled by gas spring

- Lowest height: $600\text{mm} \pm 10\text{mm}$
- Highest height: $920\text{mm} \pm 10\text{mm}$
- Trendelenburg tilting back & forth is controlled by the foot pedal on the sides of the stretcher.
- Safety working load: 2500N
- Using environment : temperature: $5^{\circ}\text{C} - 40^{\circ}\text{C}$
- The main frame of stretcher is made of galvanized steel tube EPOXY coated electrostatic spray.
- The top board BACKREST is made of X-ray material.