

Trauma Stretcher



Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Description - Back rest: $0-85^{\circ}\pm 5^{\circ}$, controlled by gas spring

- Lowest height: $600\text{mm}\pm 10\text{mm}$

- Highest height: $920\text{mm}\pm 10\text{mm}$

- Trendelenburg tilting back & forth is controlled by the foot pedal on the sides of the stretcher.

- Safety working load: 2500N

- Using environment : temperature: $5^{\circ}\text{--}40^{\circ}$

- The main frame of stretcher is made of galvanized steel tube EPOXY coated electrostatic spray.

- The top board BACKREST is made of X-ray material.