

## Hydraulic Stretcher



Hydraulic stretcher

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Epoxy coated mesh platform

- Air assisted backrest 0-45
- Folding bed sides.
- Height adjustment 460 - 850 mm
- 125 mm central locking castors.
- 1 V pole can be placed in all 4 corners
- Rotating buffer wheels in all corners
- Hydraulic height adjustment.
- No Tilting
- Size 2100 x 690 x 460 -850 mm