

Cruiser Wheelchair



- Size : 14" to 20"
- Frame : Steel / Folding
- Upholstery : Vinyl
- Armrest : Removable / Full
- Footrest : Swing-Away
- Footplate : ABS
- Castor : 8" x 2" Solid PU or 8" x 1"
- Rear Wheel : 24" x 1.3/8" Mag Solid PU
- Max. Carry Capacity : 110kgs
- Weight of Chair : 18.8kgs
- Back Height : 42cm
- Seat Depth : 40cm

Rating: Not Rated Yet

Price

[Ask a question about this product](#)