

## Manual Steel Wheelchair



- Size : 14" , 18" & 20"
- Frame : Steel / Folding
- Upholstery : Nylon
- Armrest : Removable / Desk
- Footrest : Removable / Steel hanger
- Footplate : ABS
- Castor : 8" x 2" Solid PU
- Rear wheel : 24" x 1.3/8" Spoke Solid PU
- Max. Carry Capacity : 80kgs – 90kgs
- Weight of Chair : 21kgs

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)