

## Trauma Stretcher



Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Back rest:  $0-85^{\circ}\pm 5^{\circ}$ , controlled by gas spring

- Lowest height:  $600\text{mm}\pm 10\text{mm}$
- Highest height:  $920\text{mm}\pm 10\text{mm}$
- Trendelenburg tilting back & forth is controlled by the foot pedal on the sides of the stretcher.
- Safety working load: 2500N
- Using environment : temperature:  $5^{\circ}\text{--}40^{\circ}$
- The main frame of stretcher is made of galvanized steel tube EPOXY coated electrostatic spray.
- The top board BACKREST is made of X-ray material.